



Temperature Scanner FAQ's

1. When am I supposed to get my temperature scanned?
 - a. Every Bilsteiner should have their temperature scanned immediately upon entering the plant and should be taken once a day.
2. Where do I get my temperature scanned?
 - a. There are 2 scanners set up. One in the AM Breakroom and one in OE Breakroom.
3. How do I get my temperature scanned?
 - a. Remove any hats, glasses and mask; stand in the green-taped box on the floor next to the scanner; once the monitor beeps, scan your badge and take your sticker.
4. Am I allowed to scan my temperature more than once per day?
 - a. Yes. If you'd like to re-scan before leaving, to make sure you are safe to go home to your family, you are welcome to do so.
5. What happens if I fail the screening?
 - a. If you do not pass the first screening, reapply your mask and glasses, sit in the marked chair for 5 minutes then recheck your temperature at the scanner (don't forget to remove your mask and any glasses/hats).
 - i. If you pass the second time, you can proceed to scan your badge and get your sticker
 - ii. If you fail the screening a second time, use the thermometer to take your temperature and follow the instructions printed at the thermometer station.
 1. If your temperature is below 100.4 degrees, you can reapply your mask and glasses, notify your supervisor and proceed to work
 2. If your temperature reads 100.4 or higher, reapply your mask, immediately leave the building and return home. Call your primary care provider (or TriHealth Priority Care at 513-346-3399 if you do not have a primary care provider) *and* call Nurse Teri (513-341-0659).